


Application of Patient-Derived Induced Pluripotent Stem Cells (iPSCs) to Study the Role of Neurite Inhibition and Mechanisms of Recovery in Alzheimer’s Disease

K. Clarke¹, M. Roger¹ and S. Przyborski^{1,2}

¹*Department of Biosciences, Durham University, South Road, Durham, UK, DH1 3LE.*
²*Reprocell Europe, NETPark Incubator, Sedgefield, UK, TS21 3FD.*



Introduction

Alzheimer’s disease (AD) is a complex neurodegenerative disease that leads to a clinical decline in neurological function due to a **loss of neuronal connectivity** in the brain¹⁻³. The aberrant processing and accumulation of proteins is a hallmark of AD, particularly the extracellular **deposition of β -amyloid ($A\beta$)** peptides that leads to **senile plaque** formation and **morphological abnormalities** in developing neurites⁴. An increase in the ratio of the 42 amino acid $A\beta$ species ($A\beta_{42}$), that is longer and more insoluble, compared with the 40 amino acid long $A\beta$ species ($A\beta_{40}$) is a feature of the disease^{5,6}.

This study utilises **iPSC-derived AD-model cell** lines to study the process of $A\beta$ induced neurite inhibition and the molecular signalling pathways that drive the process. $A\beta$ is known to bind to the **Nogo receptor (NgR)**⁷, which is involved in the inhibitory response that results in the inability of neurons to regenerate following **spinal cord injury (SCI)**^{8,9}. NgR is known to activate **Rho A and ROCK** following SCI^{8,9}, impacting actin dynamics and resulting in **neurite retraction**. Evidence also suggests that $A\beta$ induction of neurite inhibition involves Rho A and ROCK¹⁰ activation.

The principal aim of this study was to further elucidate the molecular signalling events involved in $A\beta$ -mediated neurite inhibition through adding **inhibitors of Rho A (ibuprofen)** and **ROCK (Y-27632)** signalling to **two distinct models of $A\beta$ -induced neurite inhibition**. The first of which involves the **exogenous** addition of $A\beta$ species to neuronal cultures, and the second involves the use of **iPSC-derived** patient and transfected **AD-model cells**.




Figure 1: Proposed molecular signalling pathways involved in $A\beta$ -mediated neurite inhibition. $A\beta$ is known to bind to the NgR involved in neurite inhibition following SCI. NgR activates Rho A and ROCK, which ultimately results in actin stabilisation. This effect upon cellular actin dynamics results in growth cone collapse and neurite retraction. $A\beta$ is thought to activate Rho A and ROCK, however it is not yet known if this is through activation of NgR.

Methodology




Fig 2: Neurosphere formation and induction of neurite outgrowth in 2D and 3D culture from human pluripotent stem cell lines. Neurospheres were formed from both embryonal carcinoma (EC) stem cell and induced pluripotent stem cell (iPSC) derived cell lines. The EC cell line, TERA2.cl.SP12 was first maintained in 2D monolayer culture (A) before being seeded in suspension and differentiated for 21 days with the synthetic retinoid EC23. This promoted the formation of mature neurospheres (C) for use in neurite outgrowth studies. Neurospheres were also formed from the iPSC-derived neuroprogenitor line, ReproNeuro, which were initially maintained as monolayer cultures (B). AggreWell™ plates were used to form spheroid structures (D). Neurospheres from both cell lines were subsequently seeded onto 2D (E) or 3D (F) ECM coated growth substrates for a further 10 days during which time neurites are formed. In 2D culture neurites radiate from the central neurosphere, whereas in 3D culture neurite penetrate the 200 μ m depth of the scaffold and are visible from the underside of the scaffold. To study the effect of $A\beta$ species on neurite outgrowth, exogenous $A\beta$ was added to cultures of EC-derived neurospheres and ReproNeuro AD model cell types were used to further study this mechanism.

ReproNeuro Patient-1 cells are neuroprogenitor cells derived from an AD patient with a mutation in presenilin 2 (PS2) and ReproNeuro Mutation are neuroprogenitor cells derived from a healthy individual transfected with a mutation in presenilin 1 (PS1). Presenilins form the basis of the γ -secretase complex involved in generation of $A\beta$ species and both cell types can be used to model the aberrant generation of $A\beta$ associated with AD. Scale Bars: A,B and D: 200 μ m C: 500 μ m.

Induction of Neurite Outgrowth from Healthy and AD-Model Cells





Fig 4: Neurite outgrowth is inhibited from Alzheimer's disease model cells in 2D culture and recovered by inhibition of Rho A and ROCK. ReproNeuro iPSC-derived wild type (A-C) and AD-model cells (D-I) were cultured in 2D with the selective ROCK inhibitor Y-27632 and inhibitor of Rho A ibuprofen, to identify any differences in their ability to form neurites. Patient derived neurospheres produce significantly less dense neurite outgrowth than any other cell type tested (J) and both patient and mutation AD-model cell lines produce significantly less dense neurite outgrowth than wild type cells. Medium supplementation with Y-27632 enhanced neurite density to a level greater than the control and reduced any differences between the cell types (K). Similarly, ibuprofen treatment also reduced any differences in neurite density between the cell types tested. Scale bars: 100 μ m, * = $p \leq 0.05$, ** = $p \leq 0.01$, *** = $p \leq 0.001$, **** = $p \leq 0.0001$

Neurite Inhibition Induced by Exogenous $A\beta$





Fig 3: Exogenous $A\beta$ species inhibit neurite outgrowth and is recovered through inhibition of Rho A and ROCK. Neurospheres derived from EC cells were cultured in 2D (A-D) and 3D (E-H) conditions with the exogenous addition of $A\beta_{40}$ to the culture medium (B,F). In addition to this treatment, neurospheres were also treated with the selective ROCK inhibitor Y-27632 (C,G) and inhibitor of Rho A, ibuprofen (D,H). This approach was used to determine if $A\beta_{40}$ supplementation had an inhibitory effect on neurite growth and if inhibition of ROCK and Rho A known to be involved in downstream NgR and $A\beta$ signalling was capable of recovering such an effect. Both quantification of neurite density from neurospheres cultured in 2D (I) and neurite penetration in 3D cultures (J) revealed that $A\beta_{40}$ significantly reduced neurite outgrowth, whilst Y-27632 treatment restored neurite outgrowth despite the presence of $A\beta_{40}$ enhancing levels of neuriteogenesis compared with the vehicle matched control. Ibuprofen had a similar effect in 2D culture enhancing neurite outgrowth to a level that surpassed the control, however ibuprofen treatment had no effect upon neurite outgrowth in 3D culture. Scale bars: 200 μ m, * = $p \leq 0.05$, ** = $p \leq 0.01$, *** = $p \leq 0.001$, **** = $p \leq 0.0001$

Conclusions

- Exogenous** addition of $A\beta_{40}$ to culture medium significantly **inhibits neurite outgrowth** in both 2D and 3D culture.
- Inhibition of Rho A and ROCK** by ibuprofen and Y-27632 respectively can **restore neurite outgrowth** in the presence of inhibitory $A\beta_{40}$.
- Alzheimer's disease **model cells** with common AD-associated **mutations** provide a **novel system** to study the role of $A\beta_{40}$ in neurite inhibition.
- Neurite density** in 2D culture is **reduced** in **AD-model cells** compared with their wild type counterpart.
- Inhibition of Rho A** by ibuprofen and **ROCK** by Y-27632 **restores neurite density** in AD-model cells.
- Rho A and ROCK** signalling potentially involved in **$A\beta$ -mediated neurite inhibition**.
- Neurite outgrowth assay** utilising AD-model cells can be used as a potential **screening tool** and has implications in the field of **personalised medicine**.




Fig 5: Inhibition of $A\beta$ induced signalling. It is not known whether $A\beta$ activates NgR itself, however, $A\beta$ is known to activate downstream Rho A and ROCK. This is thought to be the major signalling mechanism responsible for $A\beta$ -induced neurite inhibition. For this reason inhibition of Rho A and ROCK through modulators such as ibuprofen and Y-27632, is a potential method of intervention that could restore neurite outgrowth and neural connectivity in the context of AD.

References

1. La Ferla et al. Intracellular amyloid-beta in Alzheimer's disease. Nature Reviews Neuroscience. (2007) 8(7):599-609 2. Cummings et al. Alzheimer's disease drug-development pipeline: few candidates, frequent failures. Alzheimer's Research & Therapy. (2014) 6(4):37 3. Hardy et al. Pathways to Alzheimer's disease. Journal of Internal Medicine. (2014) 275(3):296-303 4. Grutzendler et al. Various dendritic abnormalities are associated with fibrillar amyloid deposits in Alzheimer's disease. Annals of the New York Academy of Sciences. (2007) 1097(1):30-45 5. Saave et al. Additional use of Aβ42/Aβ40 ratio with cerebrospinal fluid biomarkers P-tau and Aβ42 increases the level of evidence of Alzheimer's disease pathophysiological process in routine practice. J. Alzheimer's Disease. (2014) 41 6. Dumurgier et al. Cerebrospinal fluid amyloid-β 42/40 ratio in clinical setting of memory centers: a multicentric study. Alzheimer's Research and Therapy. (2015) 7(1):1-9 7. Park et al. Nogo Receptor Interacts with Brain APP and Aβ to Reduce Pathologic Changes in Alzheimer's Transgenic Mice. Current Alzheimer Research. (2007) 4(5):568-570 8. Yiu & He. Glial inhibition of CNS axon regeneration. Nature Reviews Neuroscience (2006) 7(8):617-626 9. Fawcett & Asher. The glial scar and central nervous system repair. Brain Research Bulletin (1999) 49 (6):377-391 10. Petratos et al. The beta-amyloid protein of Alzheimer's disease increases neuronal CRMP-2 phosphorylation by a Rho-GTP mechanism. Brain. (2008) 131:90-108