

Improvement of Human iPS Cell-Derived Hepatocyte Functionality Using 3D Culture System

Zachary Yu-Ching Lin¹, Kenichi Tamura¹, Rina Akahira¹, Robert R. Annand², Shunsuke Yoshida¹, and Yasuyuki Hayashi¹

¹ReproCELL Inc., Japan, KDX Shin-Yokohama 381 Bldg. 9F. 3-8-11, Shin-Yokohama, Kohoku-ku, Yokohama, Kanagawa 222-0033, Japan. ²Stemgent-A ReproCELL Group Company., 4 Hartwell Place, Lexington, MA 02421 USA




Introduction

Human primary hepatocytes are utilized for high-throughput screening in early-stage drug discovery in order to evaluate thousands of potential compounds. Yet, human primary hepatocytes have the disadvantage of a limited supply from a single donor as well as high donor-to-donor variability. To overcome these obstacles, functional human induced pluripotent stem (iPS) cell-derived hepatocytes are highly desirable, as they are available in unlimited quantities from the same donor. However, immaturity and donor-to-donor variability are common drawbacks of iPS cell-derived hepatocytes. To address hepatocyte maturation, we evaluated multiple methods using 3D cultivation for maturing iPS cell-derived hepatocytes. We compared different 3D culture systems with traditional 2D cultures by analyzing the expression levels of specific cytochrome P450 (CYP) enzymes that play an important role in drug-metabolism. We believed that 3D culture is able to provide a micro-environment that promotes maturation of human iPS cell-derived hepatocytes, potentially facilitating the creation of a human iPS cell-derived hepatocyte panel, which will enable assessment of donor-to-donor variability in iPS cell-derived hepatocyte function.


Materials and Methods

- ReproHepato type I™ kit (Cat. No. RCESDH001)**
- Cells 1 vial (8.25 million cells/vial)
 - Thawing Medium 1 bottle
 - Maintenance medium 1 bottle
 - Assay Medium 1 bottle
 - Supplements
- 3D culture**
- Puramatrix® (3D Matrix Inc.)
 - Nanoshuttle™ PL (n3D Bio Inc.)
 - Low attachment plate (sumitomo bakelite)
 - Bioreactor (Able-Biott 30 mL)
- RT-PCR**
- CYP3A4 TaqMan Gene Expression Assays (Life technology, Cat.No. Hs00604506_m1)
 - CYP1A2 TaqMan Gene Expression Assays (Life technology, Cat.No. Hs00167927_m1)
 - CYP2B6 TaqMan Gene Expression Assays (Life technology, Cat.No. Hs04183483_g1)
 - CYP2C9 TaqMan Gene Expression Assays (Life technology, Cat.No. Hs02383631_s1)
 - CYP2C19 TaqMan Gene Expression Assays (Life technology, Cat.No. Hs00426380_m1)
 - CYP2E1 TaqMan Gene Expression Assays (Life technology, Cat.No. Hs00559368_m1)
 - CYP2A6 TaqMan Gene Expression Assays (Life technology, Cat.No. Hs00868409_m1)
 - GAPDH TaqMan Gene Expression Assays (Life technology, Cat.No. Hs02758991gm1)
- CYP3A4 induction assay**
- Rifampicin (Sigma, Cat.No. R7382)
 - Omeprazole (Sigma, Cat.No. 104)
 - Sodium Butyrate NA (Sigma, Cat.No. 303410)
 - P450-Glo™ CYP3A4 Assay with Luciferin-IPA (Promega, Cat. No. V9002)


Characterization of ReproHepato




3D cultivation of ReproHepato ①




3D cultivation of ReproHepato ②




Hepatic differentiation by 3D culture



Optimization of Hepatic Differentiation: ①2D Preliminary Study



②Application of Preliminary Study to 3D Culture



Conclusion

Simplified FDA guidelines			
Category	CYP family	Basal expression	Induction
1	CYP3A4	similar to primary hepatocyte	at least 5 fold
		2D <input checked="" type="checkbox"/> 3D <input checked="" type="checkbox"/>	2D <input type="checkbox"/> 3D <input checked="" type="checkbox"/>
	CYP1A2	similar to primary hepatocyte	at least 5 fold
2	CYP2C9	similar to primary hepatocyte	
		2D <input type="checkbox"/> 3D <input checked="" type="checkbox"/>	
	CYP2C19	similar to primary hepatocyte	
3	CYP2D6	similar to primary hepatocyte	
		2D <input type="checkbox"/> 3D <input type="checkbox"/>	
	CYP2E1	similar to primary hepatocyte	
		2D <input checked="" type="checkbox"/> 3D <input type="checkbox"/>	
	CYP2A6	similar to primary hepatocyte	
		2D <input checked="" type="checkbox"/> 3D <input type="checkbox"/>	