

Calcium imaging of ReproNeuro cells

ReproNeuro RCDN002N, RCDN003P, RCDN001N
(Also discontinued items RCESDN001 & RCESDN002)

Version 2.0

Schedule

day	-1	0	3	7	14	...	21
Plate coating	●						
Medium preparation		●					
Thawing cells and plating		●					
Fresh Medium Exchange			●	●	●		●
Assay					●	●	●

Required Reagents and Supplies

Name	Catalog Number	Note
ReproNeuro Culture Medium	RCDN101	Store at 4°C
Hank's Balanced Salt Solution	ThermoFisher Corp, (GIBCO) Cat.No. 14025-076	Store at 4°C
Dimethyl Sulfoxide (DMSO)	Nacalai Japan Cat.No. 13407-45	
Fluo-8 [Ex/Em: 490nm/514nm]	ABD Bioquest Cat.No. 21081	Make 5mM in DMSO; store -70°C
35 mm non-coated dish with grid	Matsunami Glass Cat.No. D111505	see diagram on next page
0.01% Poly-L-Lysine (PLL) Solution	Sigma Cat.No. P4832	Store at 4°C

ReproNeuro Coat	Reprocell, RCDN201	Store at 4°C
-----------------	--------------------	--------------

Instruments for Observation

- Confocal laser scanning microscopes [Olympus, FLUOVIEW FV1000 IX81]
- Objective lens [Olympus, UPLSAPO 40 x oil]

Procedures

1. Plate ReproNeuro cells

- 1.1. Prepare 5 mL of PLL coating solution as described in **section 1** of the *Protocol of ReproNeuro™*. Coat the 35 mm dish (glass portion) with 425 uL volume of 0.0033% of PLL at 37°C for 2 hours.
- 1.2. Just prior to the next step, prepare the ReproNeuro Coat solution by mixing 1 mL of PBS and 30 uL of ReproNeuro Coat in a sterile 15 mL conical tube.
- 1.3. Retrieve the plate from the 37°C incubator. Remove the PLL solution and apply 425 uL of fresh PBS to the plate for a rinse. Remove and discard the rinse. Repeat the rinse one more time.
- 1.4. Add 425 uL volume of the previously prepared ReproNeuro Coat solution to the glass region of the plate. Incubate at 37°C overnight.
- 1.5. The next day, thaw and prepare ReproNeuro cells according to the manual, up to **step 4.12**. The cell concentration at this step is 2.0×10^5 cells/ mL. (See *Protocol of ReproNeuro™*).
- 1.6. Remove the ReproNeuro Coat solution from the plate and immediately plate 425 uL prepared ReproNeuro cells on the coated dish (cell density is 7.5×10^4 cells/cm²). Place the cell suspension only on the glass region (1.13 cm² area). Incubate covered dish at 37°C overnight for about 8-12 hours to allow attachment.

- 1.7. Remove 200 μ L of the culture medium and exchange for fresh Culture medium (200 μ L).
- 1.8. Repeat this exchange of half volume (200 μ L) of Culture medium at day 3, 7, 14, 21.
- 1.9. It is recommended to perform the calcium imaging assay between day 14 and 21.

2. Calcium Imaging Assay

- 2.1. The Fluo-8 working concentration is 5 μ M in Culture medium. To prepare the reagent, add 2 ml of Culture medium into a sterile 15 ml conical tube. Add 2 μ L of Fluo-8 (5 mM) reagent into the tube containing the Culture medium. The concentration is now 5 μ M and ready to apply to the cells.
- 2.2. Remove the cultured medium from to the dish and then add 425 μ L of the freshly prepared Fluo-8 (5 μ M) solution.
- 2.3. Incubate the dish at 37°C for 40 min.
- 2.4. Remove the Fluo-8 (5 μ M) solution from to the dish and then add 2 mL of HBSS to rinse the entire dish.
- 2.5. Remove the HBSS (salt solution rinse) from the dish and then again add 2 mL of HBSS.
- 2.6. Remove the HBSS (salt solution rinse) from the dish and then again add 2 mL of HBSS.


Note: You have now completed 3-washes

- 2.7. Incubate the dish at 37°C for 15 min and then observed the cells in the HBSS.

Matsunami Glass Industries, Ltd.

Website: <http://www.matsunami-glass.co.jp/life/index.html>

Figure 1. Schematic of the Matsunami Glass company 35mm glass bottom dish which is suitable for Calcium Imaging of ReproNeuro cells.



Special Thanks - Reprocell expresses gratitude to Mr. Enya and Dr. Oka at Kio University for providing guidance on derivation of this protocol.

ReproCELL Inc.

Website: <http://www.reprocell.com/en>

E-mail: info_en@reprocell.com